

DO IT AT HOME IMPRESSION KIT

We have a 99% success rate using our do it at home impression kit. You will be successful as well. If you go to the website you can read a bit more and have all the videos as well.

Link to the do it at home impression process

Go here to get all the info and all the videos

<https://www.bigearinc.com/providers/big-ear-home-impression-kit/>

Link to the full do it at home impression video

https://www.youtube.com/watch?time_continue=147&v=-cH152qA8rk

REPLACE ALL CONTENTS AND SHIP EVERYTHING BACK ASAP

1. Place person in comfortable chair where both ears can be easily accessed.
2. Examine ear canal using ear light. Make sure canal is clear of wax blockages.
3. Place oto block into ear canal. Make sure you go to the second bend.

The block should fit snugly and fully fill the canal. Make sure the string is hanging straight down out of the ear. Place block to depth of line on sight, be gentle, but go deep.

4. Open 2 part silicone solution.

Mix manually, rubbing together until pink throughout (15-20 seconds).

5. Roll between hands and insert into syringe.

6. Place bite block into mouth, short side side down.

For stereo earbuds use the bite block, for sleeping ear buds open your jaw, close your jaw, grind your teeth.

7. Empty syringe contents into ear, starting deep.

Be sure to fill entire ear with silicone.

8. Silicone will be firm in around 5 minutes.

9. Once firm, pull material from ear.

Start high on the ear and pull forward, twisting out of canal while holding oto block string.

10. Clean out all remaining silicone out of syringe.

11. Repeat procedure with other ear.

12. REPLACE ALL CONTENTS BACK INTO BOX EAR LITE, SYRINGE, COMPLETELY FILLED OUT FORM, AND IMPRESSIONS

13. PLACE SELF ADDRESSED LABEL ONTO BOX AND MAIL BACK TO BIG EAR

ANY QUESTIONS PLEASE CALL GLENN 719-271-9081



407 Golden Eagle Dr.
Blythewood,
SC 29016. USA
Ghood@bigearinc.com
(719) 271-9081
www.BigEarInc.com

